

# From multi-coil kspace to single-coil kspace

Zizhao Zhang

August 27, 2018

## Abstract

The way to combine multi-coil kspace data acquired from parallel MRI to a single kspace data for neural network training.

## 1 Method

Define a k-space data as  $K$  and this Fourier transformation  $\mathcal{F}^{-1}(K)$  as  $X = a + jb$ . When we have  $N$  multi-coil data, we can obtain multiple  $\{K_i\}$  and thereby  $\{X_i\}$ .

Define  $|X|$  as the magnitude of the complex data  $X$ , i.e. the MRI image we visualize. To obtain the single coil magnitude image  $|Y|$ , we compute

$$\begin{aligned} |Y| &= (|X_1|^2 + |X_2|^2 + \dots + |X_N|^2)^{\frac{1}{2}} \\ &= ((a_1^2 + b_1^2) + (a_1^2 + b_1^2) + \dots + (a_N^2 + b_N^2))^{\frac{1}{2}} \\ &= ((a_1^2 + a_2^2 + \dots + a_N^2) + (b_1^2 + b_2^2 + \dots + b_N^2))^{\frac{1}{2}} \\ &= (A^2 + B^2)^{\frac{1}{2}} \end{aligned} \tag{1}$$

where  $A = \sqrt{(a_1^2 + a_2^2 + \dots + a_N^2)}$  and  $B = \sqrt{(b_1^2 + b_2^2 + \dots + b_N^2)}$ .

Thus, we can treat

$$Y = A + jB \tag{2}$$

as the single coil image data obtained from multi-coil data without losing information. And we can compute the corresponding kspace data as  $\mathcal{F}(Y)$ . Note that the sign the real and imaginary parts can not be retrieved. See the visual result below.




Figure 1: **magnitude**: The standard method to generate an image from multi-coil kspace data. **method magnitude**: The proposed method to compute magnitude of  $Y$ . **method image real**: The real part of  $Y$ . **method image imaginary**: The imaginary part of  $Y$ . **method kspace magnitude**: The magnitude of  $\mathcal{F}(Y)$ . **method kspace phase**: The phase of  $\mathcal{F}(Y)$ .